

EasyAid

ACTIOMEDIC®




EasyAid


The innovative plaster dispenser

from ACTIOMEDIC®

Schnelle Hilfe in nur wenigen Schritten



Pull the plaster out towards the bottom.



Place plaster onto skin.
Pull off second half of
sealing paper.



Apply plaster and smooth
down firmly.

Discover the EasyAid plaster dispenser now!

easyaid.gramm-medical.de

Clean & tidy

- ✓ removal of the plasters out towards the bottom
- ✓ no dirt or blood drops in the dispenser or on remaining plasters



Anti-theft lock

- ✓ optional activatable anti-theft lock
- ✓ dispenser housing only opens with the matching key



The EasyAid Reminder

- ✓ implemented, smart Kanban system
- ✓ on time reminder to reorder new plaster refills



Well-protected

- ✓ removal of the sealing paper only upon use
- ✓ completely closed housing
- ✓ efficient protection against dirt, dust and splash water



3-fold protection

- ✓ offers room for up to three different plaster varieties



Easy to find


- ✓ signal color orange
- ✓ very good visibility
- ✓ high discernibility even in challenging environments



EasyAid. Easy life.


Sophisticated design meets uncomplicated handling and efficient functionality.

The EasyAid plaster dispenser is mainly defined by its effortless handling and smart design. In contrast to other plaster dispensers, the EasyAid has a completely closed housing. This protects its interior and its contents from dust, dirt and splash water - even in challenging environments. After opening the transparent protective lid,




the plasters can be pulled out straight down. This construction, along with the fact that the sealing paper is only removed when the plaster is pulled out of the dispenser, keeps the remaining plasters well-protected from dirty hands and dripping blood at all times. Its signal color makes it easy to find,

turning it into a central and visible emergency contact point. The EasyAid plaster dispenser also benefits from an uncomplicated



and quick installation - in many cases the drill holes of an old plaster dispenser can even be reused (applies to most current models). Alternatively, adhesive strips can be used to mount the dispenser. With slots for up to three plaster refills, the EasyAid offers sufficient space for different types of plasters and can therefore be individually adapted to its field of application. The EasyAid plaster dispenser can also



be used by large companies and public institutions thanks to an optionally activatable anti-theft lock. Another highlight: the integrated Kanban system. Using a reminder directly implemented in the plaster refill, you will always be informed on

time to reorder new refills.

Say goodbye to always empty plaster sets and dispensers - say hello to EasyAid!

EasyAid


The innovative plaster dispenser by ACTIOMEDIC®



EasyAid plaster dispenser, Standard I

60x **AQUATIC** plaster strips 19 x 72 mm
30x **ELASTIC** plaster strips 19 x 72 mm

Order Nr. 120.800.10000 | Qty: 1 pcs.



EasyAid plaster dispenser, Standard II

30x **ELASTIC** plaster strips 19 x 72 mm
30x **AQUATIC** plaster strips 19 x 72 mm
15x **ELASTIC** finger tip plasters 43 x 68 mm


Order Nr. 120.800.12000 | Qty: 1 pcs.



EasyAid plaster dispenser, Standard III

30x **ELASTIC** plaster strips 19 x 72 mm
24x **ELASTIC** plaster strips 120 x 20 mm
15x **ELASTIC** finger tip plasters 43 x 68 mm


Order Nr. 120.800.14000 | Qty: 1 pcs.



EasyAid plaster dispenser, KIDS I

90x **KIDS** plaster strips 19 x 72 mm

Order Nr. 120.800.30000 | Qty: 1 pcs.



EasyAid plaster dispenser, Detect I

60x **DETECT AQUATIC** plaster strips 19 x 72 mm
30x **DETECT ELASTIC** plaster strips 19 x 72 mm


Order Nr. 120.800.20000 | Qty: 1 pcs.



EasyAid plaster dispenser, Detect II


30x **DETECT AQUATIC** plaster strips 19 x 72 mm
30x **DETECT ELASTIC** plaster strips 19 x 72 mm
24x **DETECT ELASTIC** plaster strips 120 x 20 mm

Order Nr. 120.800.22000 | Qty: 1 pcs.



Plaster Refills

for the EasyAid plaster dispenser



30x DETECT AQUATIC plaster strips, 19 x 72 mm

metal detectable, water-repellent & hypoallergenic

Order Nr. 120.823.19720 | Qty: 8 pcs.

30x DETECT ELASTIC plaster strips, 19 x 72 mm

metal detectable, elastic & hypoallergenic

Order Nr. 120.824.19720 | Qty: 8 pcs.



24x DETECT AQUATIC finger plaster, 120 x 20 mm

metal detectable, elastic & hypoallergenic

Order Nr. 120.824.12020 | Qty: 8 pcs.



30x AQUATIC plaster strips, 19 x 72 mm

water-repellent, hypoallergenic & breathable

Order Nr. 120.821.19720 | Qty: 8 pcs.



30x ELASTIC plaster strips, 19 x 72 mm

elastic, hypoallergenic & air-permeable

Order Nr. 120.820.19720 | VE: 8 pcs.

15x ELASTIC finger tip plasters, 68 x 43 mm

elastic, hypoallergenic & air-permeable

Order Nr. 120.821.04086 | Qty: 8 pcs.



24x ELASTIC finger plaster, 120 x 20 mm

elastic, hypoallergenic & air-permeable

Order Nr. 120.820.12020 | Qty: 8 pcs.



30x KIDS plaster strips, 19 x 72 mm

water-repellent, hypoallergenic & breathable

Order Nr. 120.826.19720 | Qty: 8 pcs.